



AFRL-RX-TY-TP-2008-4617

PREPRINT

THE DRY AEROSOL DEPOSITION DEVICE (DADD): AN INSTRUMENT FOR DEPOSITING MICROBIAL AEROSOLS ONTO SURFACES

Brian K. Heimbuch, Kimberly R. Kinney, Robert K. Nichols
Applied Research Associates
P.O. Box 40128
Tyndall Air Force Base, FL 32403

Joseph D. Wander
Air Force Research Laboratory

DECEMBER 2008

**Distribution Statement A:
Approved for public release; distribution unlimited.**

This document was submitted for publication in the Journal of Microbial Methods.

**This material is declared a work of the U.S. Government and is not subject to
copyright protection in the United States.**

**AIRBASE TECHNOLOGIES DIVISION
MATERIALS AND MANUFACTURING DIRECTORATE
AIR FORCE RESEARCH LABORATORY
AIR FORCE MATERIEL COMMAND
139 BARNES DRIVE, SUITE 2
TYNDALL AIR FORCE BASE, FL 32403-5323**

REPORT DOCUMENTATION PAGE					<i>Form Approved OMB No. 0704-0188</i>	
<small>The public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.</small>						
PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.						
1. REPORT DATE (DD-MM-YYYY)		2. REPORT TYPE			3. DATES COVERED (From - To)	
4. TITLE AND SUBTITLE				5a. CONTRACT NUMBER		
				5b. GRANT NUMBER		
				5c. PROGRAM ELEMENT NUMBER		
6. AUTHOR(S)				5d. PROJECT NUMBER		
				5e. TASK NUMBER		
				5f. WORK UNIT NUMBER		
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES)					8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)					10. SPONSOR/MONITOR'S ACRONYM(S)	
					11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT						
13. SUPPLEMENTARY NOTES						
14. ABSTRACT						
15. SUBJECT TERMS						
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT	18. NUMBER OF PAGES	19a. NAME OF RESPONSIBLE PERSON	
a. REPORT	b. ABSTRACT	c. THIS PAGE			19b. TELEPHONE NUMBER (Include area code)	

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23

Title

The Dry Aerosol Deposition Device (DADD):
An Instrument for Depositing Microbial Aerosols onto Surfaces

Authors and affiliation

Heimbuch, B.K., Kinney, K., Nichols, B.
Applied Research Associates, 139 Barnes Drive, Suite 2,
Tyndall Air Force Base, Fla. 32403

Wander, J.D.
Air Force Research Laboratory, Airbase Technologies Division (RXQ)
139 Barnes Drive, Suite 2, Tyndall Air Force Base, Fla. 32403-5323

Author contact information

Brian Heimbuch	850-283-9724	brian.heimbuch.ctr@tyndall.af.mil
	Fax: 850-283-9707	
Kimberly Kinney	850-283-4025	kimberly.kinney.ctr@tyndall.af.mil
Bob Nichols	850-283-6061	bob.nichols.ctr@tyndall.af.mil
Joe Wander	850-283-6240	joe.wander@tyndall.af.mil

24

Abstract

25 Contamination of infrastructure and equipment with biowarfare agents has led to the
26 development of antimicrobial surfaces/coatings that are designed to “self-sterilize.”
27 Surfaces will likely be contaminated via an aerosol exposure and thus antimicrobial
28 efficacy measurements should also be performed using biological aerosols. Standard
29 tests that use microbial agents suspended in aqueous buffers may provide misleading
30 results that overestimate the performance of the surface. A settling chamber is a
31 common instrument for applying biological aerosols to surfaces. However, settling
32 chambers have some drawbacks that make them undesirable for all applications (i.e.,
33 slow loading times, large footprint, variable loading, etc.). We developed a Dry Aerosol
34 Deposition Device (DADD) that uses impaction rather than settling for loading surfaces
35 with biological aerosols. The use of impaction allows for rapid and highly reproducible
36 loading of microorganisms onto surfaces. We demonstrated that the DADD can deliver
37 both *Bacillus atrophaeus* spores and *Staphylococcus aureus* vegetative cells to glass
38 coupons at concentrations exceeding 1×10^4 CFU/cm². The coefficient of variation (CV)
39 for sample-to-sample loading within an experiment was 13.6% for spores and 6.1% for
40 *S. aureus* cells. The DADD is also a relatively simple and inexpensive device that can
41 easily be contained within a 4-foot biological safety cabinet.

42

43

Key words

44 bioaerosol, deposition, method, microorganism, particle, spore

Introduction

The methods used to measure antimicrobial efficacy of decontamination technologies are, in principle, very simple: samples are contaminated with microorganisms, a subset is exposed to a decontaminant and a control group is not exposed to the decontaminant. Microbes are extracted and then quantified using viable plating; the antimicrobial efficacy is the ratio of the two values. Many standard tests (American Association of Textile Colorists and Chemists (AATCC) 144, American Society for Testing and Materials (ASTM) E2148, E2414) use some variation of this strategy to provide the efficacy of the decontamination under a given set of conditions. However, the conditions used to measure antimicrobial efficacy may influence the susceptibility of the microorganisms to the disinfection technologies (Betty et al., 2005, Lai et al., 2004, Prugh and Calomiris, 2006). For example, organic components are known to compete with oxidative agents (Betty et al., 2005) and reduce the effectiveness of ultraviolet light (Lai et al., 2004). Properties such as temperature, relative humidity, presence of organic material, presence of a carrier, water content, etc., must be considered to ensure the laboratory test simulates field-based exposures. Only then can laboratory data be extrapolated to field-based efficacy.

A major concern during the last decade is the deployment of biological warfare agents. Although the primary effect of biological agents is causing infections, a secondary consequence is contamination of infrastructure. This was clearly demonstrated during the 2001 anthrax attacks, in which multiple buildings were closed for months to allow for decontamination (Canter, 2004). Biological warfare agents are expected to be delivered

in aerosol form; therefore, a biological aerosol is an appropriate challenge to evaluate decontamination technologies. Biological aerosols are complex entities and changes in their particle size or the presence of inert materials may influence their susceptibility to the decontaminant. Another important factor affecting biological aerosols is their low water content. Water quickly evaporates from aerosol droplets, leaving droplet nuclei. The droplet nuclei are responsible for contaminating surfaces; thus droplet nuclei must be used to evaluate decontamination technologies.

No standard test methods exist for applying biological aerosols to surfaces. The likely reason is that the application of biological aerosols to surfaces is very challenging. In liquid systems, it is quite simple to apply a standard inoculum of microorganisms to surfaces. However, quantifiably adding biological aerosols to surfaces requires special instrumentation and expertise in aerosols. The commonest device for applying aerosols to surfaces is a settling chamber (Barron et al., 2008, Brown et al., 2007, Feather and Chen, 2003, Marple and Rubow, 1983, McCready, 1986). A settling chamber provides a realistic challenge and, if controlled properly, can produce reproducible results. However, one drawback to using a settling chamber is the time required to load samples. Most protocols call for loading times ranging from 10–24 hours and, depending on the requirements for the challenge, this time may be unacceptable. Also larger particles settle more quickly than smaller particles, which may cause problems when evaluating self-decontaminating surfaces due to varying contact times. Static charge on the particles or the surfaces must also be considered when loading microbes in a settling chamber as charge may affect loading and distribution of the

microorganisms (Lai, 2006). To simplify loading of samples with microorganisms, we developed a system that utilizes impaction rather than settling: The Dry Aerosol Deposition Device (DADD), which is much smaller than a conventional settling chamber and allows for rapid and highly reproducible loading of samples. This report describes the design, operation, and validation of the DADD.

Dry Aerosol Deposition Device (DADD) - Description

The DADD (Figure 1) was designed to fit in a 4-foot biological safety cabinet (BSC) to contain fugitive aerosols. The aerosol is created using a single-jet Collison nebulizer (BGI Inc, Waltham, Mass.), which produces an airflow of ~ 2 liters per minute (LPM). The aerosol passes through a diffusion dryer (TSI Inc., Minneapolis, Minn.) that uses silica gel to dry the aerosol. The silica gel surrounds the circumference of the dryer, resulting in very little particle disruption. A relative humidity of 20% is achieved in the air exiting the drier. The aerosol then flows into a two-way valve (GC Valves, Charlotte, N.C.) that is controlled by an electronic switch with a sensitivity of 0.1 second. The valve defaults to the overflow position, which directs the aerosol to a high-efficiency particulate air (HEPA) filter. The collector is a single-stage (d_{50} 0.65 μm) cascade impactor (Tisch Environmental, Cleves, Ohio) that was modified to include a motor that rotates the samples at 20 rpm during loading (Figure 2). The impactor is connected to a vacuum pump that draws air at 1 cubic foot per minute (CFM). The collector is not directly attached to the two-way valve, which allows make-up air to enter the sampler at the junction. The gap is required to account for the difference in flow rates between aerosol generator (2 LPM) and the collector

(28.3 LPM). Since the device is operated in a BSC, the make-up air is HEPA- filtered and will not contain contaminants. The two-way valve and the rotating table are both activated by the same switch, which triggers the rotating to halt once sampling is complete.

Materials and Methods

Microorganisms: *Bacillus atrophaeus* (Bg) spores (ATCC 9372) were prepared to > 95% purity using standard protocols (Nicholson and Setlow, 1990). The spores were stored at -80 °C for long-term storage and 4 °C for short-term storage. The spores were analyzed by phase contrast microscopy to ensure they were phase bright prior to use. *Staphylococcus aureus* (ATCC 6538) was prepared by growing an overnight culture in trypticase soy broth (Becton Dickinson, Franklin Lakes, N.J.) at 37 °C/220 rpm. The cells were harvested by centrifugation (10 min at 10,000 X g), then resuspended in 1% raffinose to an OD₅₅₀ = 0.8.

Test Substrates: Glass slides were cut into 1-inch square samples, washed with ethanol, and stored in sterile containers until needed.

Aerosol exposure: Bg spores were diluted in sterile water to a concentration of 1×10^7 CFU/mL and 35 mL of the spore solution was added to the single-jet Collison nebulizer. Compressed air (20 psi) was added to the Collison nebulizer and allowed to equilibrate for 5 minutes. Glass coupons were placed in glass Petri dishes (Tisch Environmental, Cleves, Ohio), which were loaded into the collector. Care was taken not to load coupons in the center of the dish or too close to the edge as the impactor does not contain jets in these areas. The vacuum source was turned on and the DADD was

activated to initiate sampling. The sampling was carried out for various times, after which the coupons were transferred into 50-mL centrifuge tubes containing 35 mL of neutralizing lecithin buffer (1X phosphate buffered saline, 0.5% lecithin, 0.59% sodium thiosulfate). The spores were removed from the coupons by aggressive vortexing for two minutes. The solution was serially diluted and inoculated in triplicate onto TSA plates using a WASP Spiral Plater (Microbiology International, Fredrick, Md.). The TSA plates were incubated overnight at 37 °C, then enumerated using a Protocol automated colony counter (Microbiology International, Fredrick, Md.). The data were loaded into the *Prism-5* statistical analysis software package (GraphPad, 2236 Avenida de la Playa, La Jolla, California) and analyzed for variance.

The aerosol process for *S. aureus* was identical to that for Bg spores with the following exceptions: 1) the cells were diluted in 1% raffinose to an OD₅₅₀ 0.8 and 2) following loading, the samples were incubated at 0-hr and 1-hr exposures to evaluate cell death due to desiccation.

Spore distribution: Glass slides were loaded with Bg spores as previously described. The slides were observed using a stereomicroscope and an upright microscope.

Results

Bg spores loaded onto glass slides with the DADD demonstrated a linear correlation between time and loading concentration (Figure 3). The 12-minute loading is a little low but the trend is clear. This indicates that loading concentration can be tuned by adjusting loading times. Reproducibility of loading is a key factor in determining the

159 usefulness of the technique. Repetitive loading performed at a single time point (5 min)
160 over multiple days indicates the DADD is capable of reproducibly loading spores onto
161 surfaces (Figure 4). The average coefficient of variation (CV) observed for spore
162 loading within a given experiment was 13.6%. The variability is well within the
163 acceptable range of variability found in traditional antimicrobial efficacy tests.

164 Loading *S. aureus* onto surfaces is more complicated because the vegetative
165 microorganism is not as hardy as a spore and may be injured during impaction or by
166 desiccation. To prevent desiccation, *S. aureus* was aerosolized in a solution of 1%
167 raffinose. *S. aureus* co-aerosolized with 1% raffinose provides a high loading
168 consistency (Figure 5). The average CV for loading triplicate samples was only 6.1%.
169 These data demonstrate that *S. aureus* survives the impaction process. As part of this
170 experiment, a duplicate set of glass slides was loaded with *S. aureus*, incubated at
171 room temperature for 1 hour, and extracted. The data from this test demonstrate that,
172 during the incubation period, an average of 13% of the cells died due to desiccation
173 compared to control samples (Figure 6). For decontamination studies, the decrease in
174 viability would be observed in both the test and the control population and would not
175 bias the result. The important factor is that a significant majority of the cells remain
176 viable during the 1-hour incubation period that is required to evaluate the performance
177 of the decontamination technology.

178 Microscopic examination of Bg spores deposited onto glass slides revealed a concentric
179 circular distribution pattern (Figure 7). The pattern is a direct result of the configuration
180 of the impactor plate, which contains a series of jets aligned in concentric circles. The *S.*
181 *aureus* cells were distributed in the same pattern as the Bg spores (data not shown).

182

183

Discussion

184 The DADD provides a rapid, highly reproducible means for challenging surfaces with
185 aerosolized microorganisms. The primary reason for developing the DADD was to
186 challenge antimicrobial/self-decontaminating/reactive materials with dry biological
187 aerosols, as water may increase the effectiveness of these surfaces. The DADD is also
188 well suited to contaminate carriers that are subsequently used to evaluate the efficacy
189 of external decontamination technologies such as liquid disinfectants or UV light. The
190 absence of water in the loading step is a key parameter in simulating threat-
191 representative aerosols. Water rapidly evaporates from aerosolized droplets (manmade
192 and naturally occurring) leaving droplet nuclei. The droplet nuclei are what eventually
193 contacts the surface and will require disinfection. A settling chamber could provide the
194 same result; however, the settling chamber requires additional time to load samples and
195 loading variability can be difficult to control.

196 The use of impaction to collect microorganisms onto surfaces is not new. The cascade
197 impactor is a common device that has been used for decades to quantify the airborne
198 concentration of particles. However, the use of the impactor for this application is unique
199 in that no attempt is being made to measure the airborne concentration; rather, the
200 impactor is being used simply to load surfaces with a quantifiable amount of
201 microorganisms. The cascade impactor typically uses collection media that are soft in
202 nature (either agar-based media or filters) or oil-coated surfaces, which limit the amount
203 of particle bounce, allowing for accurate sampling. For this study, a hard surface, glass,
204 was used for the collection of microorganisms because the glass slide can easily be

205 observed with a light microscope. Undoubtedly, many particles simply bounce off the
206 surface and are not collected. However, no attempt was made to measure aerosol
207 concentration so particle bounce is of little concern for this exercise.

208 One potential drawback to the DADD methodology is that the spores are not uniformly
209 distributed across the entire surface. Microscopic examination of Bg spores deposited
210 onto glass slides shows a concentric circular distribution (Figure 7). The pattern is a
211 direct result of the configuration of the impactor plate, which contains a series of jets
212 aligned in concentric circles. The DADD could also be used without the rotating disc,
213 which would deposit a series of evenly distributed spots on the coupon. In either case,
214 the microbes are distributed in a standard pattern over the surface, but higher
215 concentrations exist at the loading sites. Complete dispersal over the surface would be
216 desirable and the DADD could be modified to increase distribution by varying the
217 sample location under the jets. However, homogeneous distribution is not required to
218 evaluate antimicrobial efficacy of the samples. The microorganisms are deposited as
219 single particles, so their exposure to decontaminants (self contained or external) will not
220 be affected by the loading pattern.

221 The DADD is a versatile device that allows for creation and deposition of aerosols onto
222 surfaces. The characteristics of the aerosol can be changed by altering the composition
223 of the nebulization fluid or by changing the atomizer to create larger droplets. By altering
224 the composition of inert components in the aerosolization fluid, a specific threat can be
225 approximated (*i.e.*, biowarfare release, respiratory transmission, etc.). For this initial
226 study, we did not focus on trying to mimic a given threat; instead, we focused on
227 demonstrating the reproducibility of the method.

Another important factor in the evaluation of disinfection technologies is the level of microbial agglomeration. Agglomerated microorganisms will produce a more rigorous challenge because the exterior microbes will shield the interior microbes from the decontamination agents. The DADD may provide a mechanism for loading various-sized agglomerates onto surfaces. The impactor uses plates with different-sized jets for depositing different-sized particles onto surfaces. For this study, the plate with a d_{50} of 0.65 μm was used. The d_{50} is a measure of the collection efficiency of the sampler and indicates the particle size at which the sampler has a collection efficiency of 50% (Jensen et al., 1992). The cascade impactor also has impactor plates for collecting particles ranging up to 9 μm . By using different plates, it may be possible to collect agglomerates of specific sizes. We are not aware of any other technology that can be used to load agglomerates of a specific size onto surfaces.

Summary

The DADD provides a mechanism to load coupons with a highly reproducible challenge of microorganisms. The challenge is rapid and, consequently, offers an advantage over a settling chamber. The distribution pattern of microbes onto surfaces is not completely uniform, but the DADD could be modified to increase surface distribution. The data generated in this study were based on glass coupons, but we have preliminary data indicating that the device can also be used to load such other substrates as aluminum, concrete, and fabrics. Thus the use of the DADD is limited only to what can be placed in the cascade impactor. Future studies will focus on loading microorganisms onto multiple surfaces and determining the conditions required to load specific-sized

agglomerates onto surfaces. The study of agglomerates is more difficult and the DADD may provide a unique capability for understanding how agglomerated microorganisms react with decontamination technologies.

References

- AATCC 100 *Antimicrobial Finishes on Textile Materials: Assessment of*, vol. TM 100-1999. American Association of Textile Colorists and Chemists (AATCC), Research Triangle Park, N.C. 27709.
- ASTM E 2414 -05: *Standard Test Method for Quantitative Sporicidal Three-Step Method (TSM) to Determine Sporicidal Efficacy of Liquids, Liquid Sprays, and Vapor or Gasses on Contaminated Carrier Surfaces*, vol. E 2414 - 05. American Society for Testing and Materials (ASTM) International, 100 Barr Harbor Drive, West Conshohocken, Pa. 19428-2959.
- ASTM E2180 *Standard Test Method for Determining the Activity of Incorporated Antimicrobial Agents in Polymeric or Hydrophobic Materials*, vol. E 2180-01. American Society for Testing and Materials (ASTM) International, 100 Barr Harbor Drive, West Conshohocken, Pa. 19428-2959.
- Barron, P. A., C. F. Estill, G. J. Deye, M. J. Hein, J. K. Beard, L. D. Larsen, and G. E. Dahlstrom. 2008. Development of an Aerosol System for Uniformly Depositing *Bacillus anthracis* Spore Particles on Surfaces. *Aerosol Sci. Tech.* 42:159–172.

272 Betty, R. G., J. M. Bieker, and M. D. Tucker. 2005. *Agricultural Pathogen*
273 *Decontamination Technology—Reducing the Threat of Infectious Agent Spread*, vol.
274 SANDIA REPORT, SAND2006-0182 - Unlimited Release. Sandia National
275 Laboratories, Albuquerque, N.M. 87185.

276 Brown, G. S., R. G. Betty, J. E. Brockmann, D. A. Lucero, C. A. Souza, K. S. Walsh, R.
277 M. Boucher, M. Tezak, M. C. Wilson, and T. Rudolph. 2007. Evaluation of a Wipe
278 Surface Sample Method for Collection of *Bacillus* Spores from Nonporous surfaces.
279 *Appl. Environ. Microb.* 73:706–710.

280 Canter, D. A. 2004. Remediating anthrax-contaminated sites: Learning from the past to
281 protect the future. *Chem. Health Saf.* 12:13–19.

282 Feather, G. A., and B. T. Chen. 2003. Design and use of a Settling Chamber for
283 Sampler Evaluation under Calm Conditions. *Aerosol Sci. Tech.* 37:261–270.

284 Jensen, P. A., W. F. Todd, G. N. Davis, and P. V. Scarpino. 1992. Evaluation of Eight
285 Bioaerosol Samplers Challenged with Aerosols of Free Bacteria. *Am. Ind. Hyg. Assoc.*
286 *J.* 53:660–667.

287 Lai, A. C. K. 2006. Investigation of electrostatic forces on particle deposition in a test
288 chamber. *Indoor Built Environ.* 15:179–186.

289 Lai, K. M., H. A. Burge, and M. W. First. 2004. Size and UV Germicidal Irradiation
290 Susceptibility of *Serratia marcescens* when Aerosolized from Different Suspending
291 Media. *Appl. Environ. Microb.* 70:2021–2027.

292 Marple, V. A., and K. L. Rubow. 1983. An Aerosol Chamber for Instrument Evaluation
293 and Calibration. *Am. Ind. Hyg. Assoc. J.* 44:361–367.

294 McCready, D. I. 1986. Wind Tunnel Modeling of Small Particle Deposition. *Aerosol Sci.*
295 *Tech.* 5:301–312.

296 Nicholson, W. L., and P. Setlow. 1990. Sporulation, Germination and Outgrowth, p.
297 391–450. *In* C. R. Harwood and S. M. Cutting (ed.), *Molecular Biology Methods for*
298 *Bacillus*. John Wiley and Sons.

299 Prugh, A., and J. J. Calomiris. 2006. *Inactivation of Bacillus anthracis Spores Delivered*
300 *as Liquid Suspension or Aerosol to Self-Decontaminating Fabric*. Defense Technical
301 Information Center (DTIC) AFRL-HE-WP-TP-2006-0060.

302

303

304




Figure 1: The Dry Aerosol Deposition Device






Figure 2: The DADD collector is composed of a single-stage cascade impactor that was modified to contain a spinning disc that rotates the samples during loading



317

318 Figure 3: Recovery of *Bacillus atrophaeus* spores loaded onto glass coupons using the
 319 Dry Aerosol Deposition Device

320



321

322 Figure 4: Recovery of *Bacillus atrophaeus* spores loaded onto glass coupons using the
 323 Dry Aerosol Deposition Device. Loading time was five minutes.




Figure 5: Recovery of *Staphylococcus aureus* loaded onto glass coupons using the Dry Aerosol Deposition Device. Loading time was five minutes.





Figure 6: Recovery of *Staphylococcus aureus* loaded onto glass coupons using the Dry Aerosol Deposition Device. Loading time was five minutes.



332

333

334 Figure 7: Distribution of *Bacillus atrophaeus* spores loaded onto Glass Coupons using
335 the Dry Aerosol Deposition Device.

336

337